

# Isothermal titration calorimetry: Principles and experimental design



imagination at work

# Agenda

Overview of Isothermal Titration Calorimetry

ITC experimental design

Data analysis

Troubleshooting

# What is isothermal titration calorimetry (ITC)

A direct measurement of the heat generated or absorbed when molecules interact

# Microcalorimetry offers enhanced information content

Label-free

In-solution


No molecular weight limitations

Optical clarity unimportant

Minimal or no assay development

# How Do ITCs Work?

The DP is a measured power differential between the reference and sample cells necessary to maintain their temperature difference at close to zero



# Performing an ITC experiment

Ligand in syringe

Macromolecule in sample cell


Heat of interaction is measured

Parameters measured from a single ITC experiment:


Affinity -  $K_D$

Energy (Enthalpy) -  $\Delta H$

Number of binding sites -  $n$




# ITC – Before titration




- Ligand – in syringe
- ~ Macromolecule in ITC cell

# Titration begins: First injection





- Ligand in syringe
- ~ Macromolecule in cell
- ~ Macromolecule-ligand complex



As the first injection is made, all injected ligand is bound to target macromolecule.





# Return to baseline





The signal returns to baseline before the next injection.

# Second injection





As a second injection is made, again all injected ligand becomes bound to the target.

# Second return to baseline





Signal again  
returns to baseline  
before next  
injection.

# Injections continue





As the injections continue, the target becomes saturated with ligand, so less binding occurs and the heat change starts to decrease.

# Injections continue





As the injections continue, the target becomes saturated with ligand so less binding occurs and the heat change starts to decrease.

# End of titration






When the macromolecule is saturated with ligand, no more binding occurs, and only heat of dilution is observed.

# Experimental results



# MicroCal™ ITC systems

MicroCal™ VP-ITC	MicroCal™ iTC <sub>200</sub>	MicroCal™ Auto-iTC <sub>200</sub>
<ul style="list-style-type: none"><li>• 1400 <math>\mu</math>L cell</li><li>• Manual sample loading</li><li>• Up to 5 samples/day</li></ul>	<ul style="list-style-type: none"><li>• Sensitive</li><li>• Fast</li><li>• Easy to use</li><li>• <math>K_D</math> from mM to nM</li><li>• 200 <math>\mu</math>L cell</li><li>• Upgradable to full automation</li></ul>	<ul style="list-style-type: none"><li>• Unattended operation</li><li>• Up to 75 samples/day (<i>using single injection method</i>)</li><li>• <math>K_D</math> from mM to nM</li><li>• Sample cell is 200 <math>\mu</math>L</li><li>• Easy to use</li><li>• 96-well plate format</li></ul>
 A compact, light-colored MicroCal VP-ITC instrument with a single injection syringe and a large sample cell on top.	 A compact, light-colored MicroCal iTC-200 instrument with a syringe and a sample cell, similar to the VP-ITC but with a different internal configuration.	 A larger, more complex MicroCal Auto-iTC-200 instrument with a front-loading sample cell and a control panel on the right side.



# Why ITC?

# Heat is a fundamental natural property...

A single titration can yield information on:

Overall binding affinity

Hydrogen bonds and van der Waals interactions

Hydrophobic and conformational effects

Stoichiometry

**calorimetry is a direct readout**

# Stoichiometry

Number of ligand binding sites per macromolecule

If one binding site the stoichiometry is 1

By convention a “Ligand” has one binding site

A “Macromolecule” can have more than one binding site

# Effective Binding Affinity Range

$K_D$  in mM to nM range

Weak binding – low C-value method

Tight binding - minimize injection volume and concentration or use competitive (displacement) binding procedure and fitting model

# Thermodynamics

$K_B$  – binding constant

$$K_D = 1/K_B = \frac{[L] \times [M]}{[ML]}$$

$$\Delta G = RT \ln K_D$$

$$\Delta G = \Delta H - T \Delta S$$

# Free energy change

- $\Delta G$  is change in free energy
- $\Delta G \leq 0$  for spontaneous process
- More negative  $\Delta G$ , higher affinity

# Enthalpy change

- $\Delta H$  – measure of the energy content of the bonds broken and created. The dominant contribution is from hydrogen bonds.
- Negative value indicates enthalpy change favoring the binding
- Solvents play a role

# Entropy change

- $\Delta S$  – positive for entropically driven reactions
- Hydrophobic interactions
- Solvation entropy (favorable) due to release of water
- Conformational degrees of freedom (unfavorable)




# Microcalorimetry provides a total picture of binding energetics

Overall binding affinity  $K_D$  correlates with  $IC_{50}$  or  $EC_{50}$ .  
This is directly related to  $\Delta G$ , the total free binding energy

$\Delta H$ , enthalpy is indication of changes in hydrogen and van der Waals bonding

$-T\Delta S$ , entropy is indication of changes in hydrophobic interaction and conformational changes


$N$ , stoichiometry indicates the ratio of ligand-to-macromolecule binding



$$\Delta G = \Delta H - T\Delta S$$

# Same affinity, different energetics! All three interactions have the same binding energy ( $\Delta G$ )


- A. Good hydrogen bonding with unfavorable conformational change
- B. Binding dominated by hydrophobic interaction
- C. Favorable hydrogen bonds and hydrophobic interaction



ITC results are used to get insights into mechanism of binding

# How to get good ITC data

# C Values



$$C = \frac{[M]}{K_D}$$

Example:


$$K_d = 100\text{nM}$$

$$[M] = 100\text{nM}, \quad C=1$$

$$[M] = 5\mu\text{M}, \quad C=50$$

$$[M]:[L] - 1:10 \text{ for } n=1$$

# C Values



# C Values in ITC

$$C = \{[M]_{\text{tot}} / K_D\} * N$$

C = 20-100 very good


C = 10-500 good

C = 1-5 and 500-1000 OK



C = < 1 and > 1000 not wanted

# ITC experimental design

## High C



## Low C



# ITC experimental design

$K_D$ (Biacore) $\mu\text{M}$	[Protein] $\mu\text{M}$	[Compound] $\mu\text{M}$	[Protein] / $K_D$
<0.5	10	100	>20
0.5-2	30	300	15-60
2-10	50	500	5-25
10-100	30	$40 \cdot K_D$	0.3-3
>100	30	$20 \cdot K_D$	<0.3

} Fixed stoichiometry





# Low C Experiments

Extending the applications of ITC

– E.g. Fragment Based Drug Discovery


Saving Protein

# C Values-Low



Note different scales

# C Values-Low



Increase ligand concentration but not the 'valuable' protein


# Competition Experiments

Extend the affinity range that ITC can be used

- Submillimolar ( $10^{-2}$ ) to picomolar ( $10^{-12}$ )

# Competition Experiments

High C Experiment  
Poor affinity estimates



# Competitive (displacement) ITC

Used to extend range of  $K_B$  determined by ITC



Tight binding ligand A and weak binding ligand B bind to same site on macromolecule




1<sup>st</sup> experiment: Ligand B titrated into macromolecule.  
Determine  $K_B$  and  $\Delta H$



2<sup>nd</sup> experiment: Macromolecule + ligand B in cell, titrated with Ligand A. Ligand A displaces ligand B.

Use displacement model for data analysis

# Tight Binders








-  Protein
-  Tight Ligand
-  Competitor



Software  
can 'pull out' the  $K_D$   
of the tight one

# Weak Binders



Software can pull the Weak  $K_D$  out



# Competition Experimental Design


$$C = [\text{cell}]/K_{D,\text{app}}$$

$$K_{D,\text{app}} = (1/K_{D,S}) / (1 + 1/K_{D,W}[W])$$

Where  $K_{D,S}$  and  $K_{D,W}$  are the affinity of the strong and weak binders respectively and  $W$  is the concentration of the weak binder


# Displacement ITC – HIV-1 Protease-Inhibitor Binding

Amprenavir




Unable to determine KB

Acetyl pepstatin



Amprenavir + acetyl pepstatin





KB of  $3.1 \times 10^{10} \text{ M}^{-1}$

*Ohtaka, et al, Protein Sci. 11, 1908-1916 (2002)*

# ITC practical considerations

# ITC practical considerations



# Good Experimental Design

## Use the correct 'C' value

$$C = [M]_{\text{tot}} \cdot n / K_D$$

C = 20-100 very good

C = 10-500 good

C = 1-5 and 500-1000 OK

C < 1 and > 1000 not wanted

# Ligand concentration


$[L] = 10 \text{ to } 20 \times [M]$

[L] - Minimum  $50 \mu\text{M}$

May need to be adjusted based on experiment

At end of ITC, for N of 1, final  $[L]/[M]$  ratio should be 2 to 4 to ensure saturation of all M binding sites

# iTC<sub>200</sub> Experimental Design Tab



# With little prior knowledge

## Good Starting Conditions

100  $\mu\text{M}$  Ligand in the Syringe and

10  $\mu\text{M}$  Macromolecule in the cell

16 x 2.5  $\mu\text{l}$  injections

Detect  $K_D$ s of 10  $\mu\text{M}$  to 10 nM

Ideal for  $K_D$ s of 2  $\mu\text{M}$  100 nM



# Sample preparation

Use dialysis or buffer exchange column

Check calibration pipettes-by weight

Retain the dialysate/exchanged buffer

Adopt and stick to a reproducible protocol

# ITC – Choice of Buffer

What buffer(s) are protein and ligand stable in? What pH?

Buffers used for other binding studies


Solubility

Requirements of additives for binding, solubility or stability

- Salt
- Detergent
- Reducing agent
- DMSO
- Other

# Choice of buffer

Buffers have ionization enthalpies:



Use buffers with  $\Delta H_{\text{ion}} \sim 0$

Including; phosphate, acetate, formate, citrate, sulfate, cacodylate, glycine

Quaternary amines (e.g. Tris) have high  $\Delta H_{\text{ion}}$

# Choice of buffer

## Avoid DTT

Unstable and undergoes oxidation

High background heat

## Use b-mercaptoethanol & TCEP

TCEP is not stable in phosphate buffer

Use conditions in which your protein is 'happy'

# Sample preparation: small molecule ligand

If ligand is too small to dialyze, be sure material is desalted prior to final preparation

Use final dialysis buffer of macromolecule to dissolve ligand

Match pH

# ITC – “Reverse Titration”

Can also have ligand in cell – lower concentration requirement.

Have macromolecule in syringe at appropriate (higher) concentration, Need to be sure can have protein at the appropriate concentration.

At saturation: no more free ligand in ITC cell

If other than 1:1 binding, use “ligand in cell” option for curve fitting. Different model used.

# iTC<sub>200</sub> Experimental design

## Injection volume and duration

- 0.5 to 3  $\mu\text{l}$  ((range 0.1-38  $\mu\text{l}$ )
- Injection rate is 0.5  $\mu\text{l}/\text{sec}$
- Spacing
- Typical 120 secs-may want to extend to 180 seconds or more if using no feedback with large heats
- Be sure baseline returns before next injection
- If this does not occur, increase time between injections

## Filter period

- 5 sec or less recommended
- Can be increased for slow reactions

# iTC<sub>200</sub> Experimental design

Typical temperature: 25 °C

Number of injections: 12-18

Reference power

- Instrument baseline
- Set 5  $\mu\text{cal}/\text{sec}$
- If very exothermic, increase setting so  $\Delta P$  does not go below zero

Initial delay

- 60 sec minimum
- Establish baseline before 1<sup>st</sup> injection



# iTC<sub>200</sub> Experimental design

## Feedback mode

- High feedback for most ITC experiments
- Low or No feedback will give better S/N but will take longer time (increase time between injections) – normally used when working with small heats

## Stirring speed

- 1000 rpm (1500 for SIM)

# Auto- iTC<sub>200</sub>

Macromolecule/protein (for cell):

Need 400  $\mu\text{l}$  in 96 well plate

Ligand (for syringe):

Need 120  $\mu\text{l}$  in 96 well plate

Set appropriate scripts for cell/syringe filling  
and cleaning

# ITC – Enthalpy Changes

$\Delta H_{\text{observed}}$  by ITC is total of :

$\Delta H_{\text{binding}}$

$\Delta H_{\text{ionization}}$


$\Delta H_{\text{conformation}}$

Any non-specific effects (buffer mismatch, pH mismatch, heat of dilution, heat of ligand dissociation)

**Need to account for these effects by appropriate controls and experimental conditions**

# Controls

## Injection of syringe material into buffer-



Peaks should be similar in magnitude to those at the end of the actual titration experiment and constant

# Controls

If they are constant-

subtract average peak size from the experimental data (Using the 'Math' command)

# Controls

If they are not constant

Buffer Mismatch –e.g. solvent component missing

Dissociation of the syringe material upon dilution into buffer- impact on apparent  $K_D$


Genuine change in heat of dilution with concentration.

# Buffer Mismatch

2% DMSO into  
2% DMSO

2% DMSO into  
1.95% DMSO


2% DMSO into  
0% DMSO



# Dissociation of Syringe Material into Buffer

Calorimetric dilution data showing the effects of different ligands on dilution of insulin

Ref: Lovatt M, Cooper A and Camillerri P (1996)  
*Eur. Biophys. J.* 24:354-357





# Changing Heats of Dilution

The control experiment should be fit to a straight line – (or geometric function (if required- Caution!))

This best fit line can then be subtracted from titration experiment before fitting to the appropriate model.

# For Each ITC Experiment

Start with clean cell and syringe


Prepare macromolecule and ligand in **matched buffer**

Perform control titration(s) to establish heat of dilution

Set appropriate scan parameters to generate full binding isotherm


# Troubleshooting

# Bent Syringe



# Syringe Height

iTC<sub>200</sub> - Syringe Holding Nut Loose




# Baseline Position/Drift


Baseline position is the first diagnostic for data quality-information on

- Cell cleanliness
- 'Sticky' proteins
- Air Bubbles
- Time between injections

# Bubbles




# Not long enough between injections







# Not long enough between injections




# 'Sticky' Proteins or Cleanliness



# 'Sticky' Proteins or Cleanliness



# 'Sticky' Proteins or Cleanliness



# Non sigmoidal binding isotherm


No Binding

No Heat

Buffer mismatch

More than one binding event

# Buffer Mismatch-No Dialysis




# Effect of DMSO Mismatch


2% DMSO into  
2% DMSO-same stock

2% DMSO into  
2% DMSO-separate stocks

2% DMSO into  
0% DMSO




# Ligand preparation from DMSO





# Match DMSO in the protein solution




# The Cure

Dialyze or use desalting column

Check for additives that are not in both cell and syringe- Also-ask what was sample purified from e.g. was protein lyophilized in buffer and not dialyzed

Check pH of final solutions-should differ by less than 0.1 pH units. This issue is common when working with high concentrations of ligand-e.g. 500  $\mu\text{M}$  and above-weak binding

# ITC: low heat



Heats for experiment same as control


# The Cure

Change experimental temperature by at least  
10 °C


AND/OR

Increase sample concentration

# Unexpected Stoichiometry



# Instrument reference power too low



Oscillating signal:  
Power below 0



# Insoluble Ligands

'One site' interactions are symmetrical and as such the ligand can be put in the cell and the protein in the syringe

# Cure-Reverse the Titration

29.5 $\mu$ M Protein titrated with  
1.1mM Compound

11.5 $\mu$ M Compound titrated with  
179 $\mu$ M Protein



<i>Parameter</i>	<i>Ligand in syringe</i>	<i>Ligand in cell</i>
n	0.99	0.97
$K_d$	104 nM	105 nM
$\Delta G^\circ$	-9.4 kcal/mol	-9.4 kcal/mol
$\Delta H^\circ$	-11.9 kcal/mol	-12.4 kcal/mol
$T\Delta S^\circ$	-2.5 kcal/mol	-3.0 kcal/mol



imagination at work Molar Ratio



# Stoichiometry

“N” is the average number of binding sites per mole of protein in your solution, assuming:

- that all binding sites are identical and independent
- that you have pure protein (and ligand)
- that you have given the correct protein and ligand concentrations
- that all your protein is active

# ITC Maintenance

Cell cleaning:

Water/buffer

20 % Contrad 70

Keep water in cell when not in use

Reference cell:

Replace water once a week

# Thank you!

Biacore and MicroCal are trademarks of GE Healthcare companies.

GE, imagination at work and GE monogram are trademarks of General Electric Company.

All third party trademarks are the property of their respective owners.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

© 2010 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden.

*All information and material contained herein are subject to copyrights owned by MicroCal, LLC, a GE Healthcare company. Any reproduction, retransmission, republication, or other use of all or part of this presentation is expressly prohibited, unless prior written permission has been granted by MicroCal, LLC. MicroCal, LLC makes no warranty, representation or guaranty as to the content, sequence, accuracy, timeliness or completeness of the information or that the information may be relied upon for any reason.*



imagination at work